

The dry-care uv - Kids edition is a good choice

Your hearing aids will be perfectly cared by comfort drying, the life span increased and repair costs reduced significantly. It is value-preserving for modern In-Ear Monitors.

Easy handling by touch on a single button, everything is handled by processor fully automatic.

In addition to drying UV-C light provides clinical cleanliness. 99,9% of bacteria, fungi and viruses will be eliminated and health risks minimized.

Due to 5V USB technology, the device can be used in hotels or cars at home or travel, plugged in any PC, TV or audio system. A standard USB 2.0 port for this is quite sufficient.

The dry-care uv - Kids edition is energy efficient, affordable in operation and safe to use.

Effective and gentle care guaranties long life hearing experience.

dry-care uv

Kids edition



Perfect drying
and care

Hygienic cleansing
by UV-C light

Mobility by
USB-standard

Dries and cares your hearing aids


Moisture and ear wax can impair the function of hearing aids and cause failures.

The dry-care uv - Kids edition offers optimal care. Comfort dry at 50°C extends life time significantly.

Ultraviolet light (253,7nm UV-C) eliminates 99,9% of bacteria, fungi and viruses.

This clinical cleanliness prevents ear infections risks and protects both your hearing aids as well as your health.

The small, portable unit offers more freedom in application.



Perfect drying by convective flow of warm air



Hygienic cleansing by UV-C light for your health



Health improving blue light



Easy handling by simple touch on a single button



Mobility by 5V USB standard



Automatic room- and temperature modulation



Automatic turn-off after comfort drying



All features controlled by microprocessor



Highest product quality with 2 year warranty



Energy-efficient, gentle and absolutely safe