

The new DRY-SPACE UV is a good choice

Thanks to the gentle, convenient drying process over 5 hours, your CI hearing implants or BAHA bone conduction implants are perfectly cared, their service life is increased and the risk of failure and repair costs are significantly reduced. It is indispensable and preserves the value of modern hearing implants.


Simple handling due to automatic one-button operation, everything else runs fully automatically controlled by a microprocessor. In addition to drying, UV-C light and blue light ensure hygienic cleanliness. Germs, bacteria and fungi are efficiently reduced by 99.9%, thus minimizing health risks.

Whether at home or on the move, the device works on a 5V-USB basis and can be used in hotels, cars, PCs, TVs or audio devices. A standard USB 2.0 port is sufficient for this. The universal USB power supply with exchangeable adapters can be used worldwide.

The DRY-SPACE UV is energy-saving, cost-effective in operation and offers safety in application. Effective and absolutely gentle care is the best guarantee for a long life with hearing implants.

DRY-SPACE UV

Electronic drying station
for cochlear implants.



Your specialist

Perfect drying and care
for bi-aural systems

Hygienic cleansing
by UV-C light

Mobility by
Universal-USB-standard



Dries and cares your Cochlear devices

Moisture, body secretions and cerumen can severely impair the function of cochlear implants and cause failures.

The DRY-SPACE UV offers optimal care for two cochlear devices. Gentle drying at 45°C (113°F) significantly prolongs the service life and reduces repairs.

Ultraviolet light (253.7nm UV-C) eliminates 99.9% of bacteria, fungi and viruses.

This clinical cleanliness prevents ear infections risks and protects both your cochlear implants as well as your health.



Perfect drying by convective flow of warm air



Hygienic cleansing by UV-C light for your health



Health improving blue light



Easy handling by simple touch on a single button



Mobility by 5V USB standard



Automatic room- and temperature modulation



Automatic turn-off after comfort drying



All features controlled by microprocessor



Highest product quality with 2 year warranty



Energy-efficient, gentle and absolutely safe