

POWERFUL VACUUM

Once thoroughly dried, the vac-clean uv 2.1® is used with its enhanced vacuum power of over -65 kPa. To do this, the attachment nozzle optimized for the type of cleaning is attached to the base head. These correspond to the LUER system, a standardized connection system for infusions and injections in the medical field. The attachment nozzles can be easily guided into openings, crevices, and sockets. The base head collects all dirt without reducing the suction power. The blue LED light (450 nm) illuminates the area to be cleaned and also has an antimicrobial hygiene function. The base head and attachment nozzles can be easily changed and cleaned in an ultrasonic bath (e.g. u-sonic 3).

The new vac-clean uv 2.1® is a medical device and meets the high requirements for medical devices (Class I) of the Medical Device Regulation (EU) 2017/745.

vac-clean uv 2.1®

Rechargeable mini vacuum cleaner
for hearing aids and hearables.




Your specialist

- ✓ Powerful vacuum pump
- ✓ LUER system standard for different applications
- ✓ Medical device according to MDR (EU)

VERSATILE

The new rechargeable mini vacuum cleaner vac-clean uv 2.1® is suitable for all hearing aids (BTE, ITE, RIC, cochlear), all charging stations, all OTC hearables, all true wireless stereo (TWS) devices, and a wide range of electronic products such as cell phones, remote controls, and much more. These products suffer from adhering dirt and moisture, which can lead to malfunctions. Continuous cleaning and care can help prevent this.



Effective suction by vacuum pump



Sonic-Brush for perfect cleansing



Blue LED light for illumination



Easy handling by only one sensor button



Rechargeable by LiPo battery



Mobility by USB-C standard



Standard LUER system adapter



All features controlled by microprocessor



Highest product quality with 2 year warranty



Energy-efficient, gentle and absolutely safe